

Non-woven fabric

Melt blown cloth

Non-woven fabric



Bacterial filtration Effectiveness  $\geq 95\%$

Non-medical  
50 PCS

GB/T 32610-2016  
EN 149:2001+A1:2009

## Disposable Protective Mask



### Instructions



1 Keep the front outside, and hang straps over ears.



2 Close to face, and stretch the mask until it fully covers the nose, mouth and chin.



3 Adjust the nose clip and press it, and keep it closely to nose and face.



4 If damaged, smelly or damp, please replace it with a new one immediately.



**España**

MASCARILLAS DESECHABLES DE 3 CAPAS



**France**

MASQUES JETABLES DE 3 COUCHES



**English**

DISPOSABLE MASKS OF 3 LAYERS



**Italiaio**

MASCHERE IGIENICHE MONOUSO A 3 STRATI

Product name: Disposable mask

Manufacturer: Anhui Zhongnan Air Defence Works Protective Co., Ltd

Address: Qianshan comprehensive economic development zone, anhui province

Packing specification: 50 pcs/ box

Shelf life: 2 years

Production date: see inner packing

**CHARACTERISTIC:**

3-layer filter not washable  
High filtration efficiency  
Integrated nasal bridge.  
Elastic bands without latex.  
Composition percentage:  
80% Non-woven fabric  
20% Meltblown fabric  
Measurements: 175\*95mm (±0.5 mm)  
Color: Blue and White...

**CARACTERÍSTICA:**

Filtro de 3 capas no lavable  
Eficiencia de filtración de High  
Puente nasal integrado.  
Bandas elásticas sin látex  
Porcentaje composición:  
80% tela no tejida  
20% tejido Meltblown  
Medidas: 175 x 95mm (±0.5mm)  
Color: Azul y Blanco...

**CARACTÉRISTIQUES:**

Filtre à 3 couches non lavable  
Efficacité de filtration élevée  
Barrette nasale intégrée  
Bandes élastiques sans latex  
Pourcentage de composition:  
80% de tissu non tissé  
20% tissu soufflé par fusion  
Dimension: 175\* 95 mm (±0.5 mm)  
Couleur: Bleu et blanc...

**CARATTERISTICA:**

Filtro a 3 strati non lavabile  
Alta efficienza di filtrazione  
Ponte nasale integrato  
Fasce elastiche senza lattice  
Percentuale di composizione:  
80% tessuto non tessuto  
20% tessuto meltblown  
Misure: 175x95mm (±0,5 mm)  
Colore: Blu e Bianco



Anti-PM2.5



Anti-bacterial



Anti-fog



Filter droplets



6 1985452 1157509

Made in china