

**M** Suitable for children aged 4-12

THREE LAYERS OF PROTECTION

# Kids Mask



BFE95 MELTBLOWN  
FILTRATION



SAFETY AND  
PROTECTION



BREATHABLE AND  
DUSTPROOF



SKIN-FRIENDLY  
AND SOFT



30PCS

## KIDS MASK

Name: Kids Mask

Packing: 30 pcs/ box

Valid Period: 2 years

Specification: 145mm \* 94mm

Usage Time: 4 hours

Production Date: See Qualified Certificate

Main Raw Materials: Nonwoven fabric, meltblown fabric, nose clip, elastic earloop.

Storage Conditions: Ambient temperature: -10℃ - 40℃, Relative humidity: ≤85%

Executive Standard: GB/T 32610-2016

Company name: Zhejiang Wanshu Environmental Sci-Tech Co., Ltd

ADD: Building 17, Xinya paper and plastic packing center, Longgang City, Wenzhou city, Zhejiang province (2nd Floor of Huahao nonwoven Co., Ltd)

## KIDS MASK

### USAGE METHOD



MAKE THE COLORED  
SIDE OUTWARDS AND  
THE NOSE BRIDGE  
STRIP UPWARDS.



PUT ON THE EARLOOP,  
MAKE THE MASK COVERS  
THE NOSE AND THE  
MOUTH COMPLETELY.



ADJUST THE NOSE  
BRIDGE STRIP, FIT THE  
MASK SNUGLY TO THE  
NOSE AND CHEEKS.