


Medizinische Einwegmaske

Anleitung:



ziehen Sie die gefaltete Maske auseinander.



befestigen Sie die Ohrschlaufen hinter dem Ohr.



streichen Sie den Nasenbügel über die Kontur der Nase für optimale Anpassung.

Notice:

Material: Vliesstoff + schmelzgeblasener Stoff (Melt-blown) + Vliesstoff

✓ 5 Masken/Packung

✓ 3-lagige medizinische OP-Maske, schützt vor Tröpfcheninfektion

✓ Größe: 17 x 9,5 cm

✓ Standard: EN14683 Typ II



Xiaotao Xiangrui Protective Equipment Co., Ltd., Group 3, Xinhe Village, Changtangkou Town Xiantao Hubei Province



2021.08.18 2 Years made in China



202108

Disposable Protective Mask (Medical)
17.5*9.5cm Typ II
Medical Device Directive (93/42/EEC)



WeMoove e.U.
Salzbergstr.82, 6067 Absam
Österreich www.physioshop24.net



EN14683-2019 Made in China

合格证 CERTIFICATE

品名: Product Name:	一次性防护口罩 Disposable protective mask
型号: Model:	平面挂耳型 Plane lug type
主要材料: Main materials:	内层35%无纺布中层30%熔喷布外层35%无纺布 Inner layer 35% non-woven fabric, middle layer 30% Melt blown fabric, outer layer 35% non-woven fabric
规格: Specifications:	17.5*9.5cm
检验员: Examination clerk:	01
执行标准: Executive standard:	T/CTCA7-2019
生产日期: Date of manufacture:	2021.09.1
有效期: Term of validity:	2年 2 years
贮存条件: Storage conditions:	阴凉通风温度-20℃-30摄氏度, 相对湿度不大于80% Cool ventilation temperature -20°C -30°C, relative humidity no more than 80%
生产厂商: Manufacturer:	安徽省桐城市奥维防护用品厂 Anhui Tongcheng Awei protective articles factory
地址: Address:	安徽省桐城市青草镇青草工业园 Qingcao Industrial Park, Qingcao Town, Tongcheng City, Anhui Province

MADE IN CHINA